VetExpert Rapid C.Brucella Ab

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ONE STEP Canine Brucella Antibody RAPID TEST KIT

VetExpert Rapid C. Brucella Ab Test

• Explanation of the test

Canine Brucellosis is an infectious disease caused by the bacteria of the genus *Brucella canis*. The VetExpert Rapid C.Brucella Ab Test is a chromatographic immunoassay for the qualitative detection of *Brucella canis* antibody in whole blood, plasma, or serum. The VetExpert Rapid C.Brucella Ab Test has a letter of T and C as "Test Line" and "Control Line" on the surface of the kit. Both the "Test Line" and "Control Line" in result window are not visible before applying any samples. The "Control Lines" is used for procedural control. Control line should always appear if this procedure is performed properly and the test reagents of control line works. A purple "Test Line" will be visible in the result window if there are enough *Brucella canis* antibodies in the specimen.

The specially selected *Brucella canis* antigen is used in test and as both capture and detector materials. These enable the VetExpert Rapid C.Brucella Ab Test to identify to *Brucella canis* antibodies in specimens, with a high degree of accuracy.

• Materials provided (2,5,10 test/kits)

VetExpert Rapid C.Brucella Ab Test contains following items to perform the assay.

- 5,10 VetExpert Rapid C.Brucella Ab Tests
- Assay Diluent bottles
- 5,10 capillary tubes for 10µl
- Instruction for use

◆ A score line for volume of 10 µl.

• Precautions

- For in-vitro diagnostic use only.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
- Avoid splashing or aerosol formation.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose all specimens, reaction kits and potentially contaminated materials, in a biohazard container.
- Do not use the test kit if the pouch is damaged or the seal is broken.

• Storage and Stability

The VetExpert Rapid C.Brucella Ab Test should be stored at room temperature. The test kit is sensitive to humidity, and as well as to heat. Perform the test immediately after removing the test device from the foil pouch. Do not use it beyond the expiration.


• Specimen Collection and Preparation

- [Whole blood]** Collect the whole blood by using the suitable anti-coagulant. Use the whole blood within 1 day after collection. Do not use the hemolysis blood.
- [Serum or Plasma]** Centrifuge whole blood to get plasma or serum specimen.
- If specimens are not immediately tested they should be refrigerated at 2 ~ 8°C. For keeping specimens more than three days with fresh, them at - 20°C or below(serum, plasma). They should be brought to room temperature before use.
- Specimens containing precipitate may yield inconsistent test results. Such specimens must be clarified before assaying.

• Procedure of the test

- Remove the test kit from the foil pouch, and place it on the flat, dry surface.
- Take 10 µl of serum, plasma, or whole blood to the dark score line of a capillary tube.
- Slowly add 10 µl of serum, plasma, or whole blood to the sample well with capillary tube with a score line for volume of 10µl and then add 2 drops (approximately 60 µl) with bottle containing diluent buffer. If the migration is not appeared within 1 minute, add one more buffer to sample well.
- As the test result, you can see the purple band in the result window of the kit.
- Interpret test results at 20 minutes.
- Please do not read after 20 minutes.

[Test Procedure]

• Interpretation of the test

- A color band which would be appeared in the left section of the result window shows that the test is works properly, and it called is the Control line (C).
- The right section of the result window called test band indicates the test results. If another color band appears on the test band. (T).

1) Negative result

The presence of only one purple color band within the result window indicates a negative result.


2) Positive result

The presence of two color bands ("T" band and "C" band) within the result window, no matter which band appears first, means positive.

3) Invalid result

If the purple color band is not visible within the result window after performing the test, the result is considered invalid. The directions may not have been followed correctly or the test may have deteriorated. It is recommended to take the test again with same specimen.

• Limitation of the test

- The VetExpert Rapid C. Brucella Ab Test will only indicate the antibody presence against *Brucella canis* in the specimen.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the veterinarian.
- If the test result is negative and clinical symptom is persist, additional testing using other clinical methods is recommended. A negative result does not preclude the possibility of canine brucellosis.

• Expected Values

The VetExpert Rapid C.Brucella Ab Test has been compared with 2- mercapto-ethanol Rapid Slide Agglutination Test. The overall accuracy is greater or equal to 90.0%.

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