VetExpert Rapid BoviD-4 Ag

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One Step BoviD-4 VetExpert Test (Research Use Only)

VetExpert Rapid BoviD-4 Ag Test

● Explanation of the test

The VetExpert Rapid BoviD-4 Ag Test is a solid phase immunochromatographic assay for the rapid, quantitative detection of Cryptosporidium antigen, Rotavirus antigen, Coronavirus antigen and E.coli k99 antigen in calf feces. The VetExpert Rapid BoviD-4 Ag Test has the letters “T” and “C” as the Test line and Control line on the surface of the device. Both the test line and control line in the result window are not visible before applying any samples. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working. A purple test line will be visible in the result window if there is enough antigen of each target in the specimen.

● Materials Provided (5 Tests/Kit)

- 1) 5/10 VetExpert Rapid BoviD-4 Ag Tests
- 2) 5/10 assay diluents
- 3) 5/10 disposable droppers
- 4) 5/10 sample collection swabs
- 5) Instruction for use

● Precautions

- 1) For best results, strict adherence to these instructions is required.
- 2) All specimens should be handled as being potentially infectious.
- 3) Do not open or remove test kits from their individually sealed pouches until immediately before their use.
- 4) Do not reuse test kit.
- 5) All reagents must be at room temperature before running the assay.
- 6) Do not use reagents beyond the stated expiration date marked on the package label.
- 7) The components in this kit have been quality control tested as a standard batch unit. Do not mix components from different lot numbers.
- 8) The assay diluent contains low concentration of sodium azide as a preservative. Sodium azide is toxic and should be handled carefully to avoid ingestion and skin contact.

● Storage and Stability

The kits can be stored at room temperature (2-30°C) or refrigerated. The test kit is stable through the expiration date marked on the package label. DO NOT FREEZE. Do not store the test kits in direct sunlight.

● Specimen Collection and Preparation


- 1) The test should be performed with the calf feces.
- 2) After collecting the specimen using swab, the specimen should be immediately extracted and tested.
- 3) If specimens are not immediately tested, they should be refrigerated at 2-8°C. For storage not less than 48 hours, freeze the specimen at -20°C or below.

● Procedure of the Test

- 1) Collect the samples from feces using the sample collection swab pre-wetted with saline solution (not provided).
- 2) Insert the sampling swab into the assay diluents tube.
- 3) Mix the sampling swab until the sample has been dissolved into the assay diluents.
- 4) Leave the tube until the large particles have settled down to the bottom of the tube. (Approximately 1 minute)
- 5) Remove the test device from the foil pouch, and place it on a flat and dry surface.
- 6) Using a disposable dropper (provided), take an aliquot from the extracted and mixed sample in the tube.

- 7) Add 4(four) drops into the each sample hole with the disposable dropper. The mixed assay diluents should be added exactly, slowly drop by drop.
- 8) As the test begins to work, you will see purple color move across the result window in the center of the test device. If the migration has not appeared after 1 minute, add one more drop of the mixed sample to the sample well.
- 9) Interpret test results at 10 minutes. Do not decide after 20 minutes.

[Test Procedure]

A color band will appear in the left section of the result window to show that the test is working properly. This band is the control band. The right section of the result window indicates the test results. If another color band appears in the right section of the result window. This band is the test band.

1) Negative result

1. Cryptosporidium antigen.
Only one band (“C”) within the result window indicates negative against Cryptosporidium antigen.


2. Rota virus antigen
Only one band (“C”) within the result window indicates negative against Rota virus antigen.

3. Corona virus antigen
Only one band (“C”) within the result window indicates negative against Corona virus antigen.

4. Escherichia coli k99 antigen
Only one band (“C”) within the result window indicates negative against Escherichia coli k99 antigen.

2) Positive result

1. Cryptosporidium antigen.
Two color bands (“T” and “C”) within the result window of Crypto Ag indicate the presence of Cryptosporidium antigen.

2. Rota virus antigen
Two color bands (“T” and “C”) within the result window of Rota Ag indicate the presence of Rota virus antigen.

3. Corona virus antigen
Two color bands (“T” and “C”) within the result window of Corona Ag indicate the presence of Corona virus antigen.


4. Escherichia coli k99 antigen
Two color bands (“T” and “C”) within the result window of E.coli k99 indicate the presence of Escherichia coli k99 antigen.

3) Invalid Result

If the purple color band is not visible within the result window after performing the test, the result is considered invalid. The directions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be re-tested.

● Limitation of the Test

- 1) This Test is for in-vitro diagnostic use only.
- 2) This test detects the presence of Cryptosporidium antigen, Rota virus antigen, Coronavirus antigen and E.coli k99 antigen in the specimen and should not be used as the sole criterion for the diagnosis of them.
- 3) As with all diagnostic tests, all results must be considered with other clinical information available to the veterinarian.
- 4) For more accuracy of immune status, additional follow-up testing using other laboratory methods is recommended.

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