### VetExpert Rapid Test FCoV Ag/Ab

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### VetExpert Rapid FCoV Ag/Ab Test Kit

- Principles


The VetExpert Rapid FCoV Ag/Ab Test Kit is a chromatographic immunoassay for the qualitative detection of Feline Corona virus (FCoV) antigens in feline feces and antibodies in feline whole blood, serum, or plasma. The VetExpert Rapid FCoV Ag/Ab Test Kit has two letters which are test(T) line and control(C) line on the surface of device. The control line is a reference line which provides the test is performed properly. The test line will appear every time when the test has performed. If FCoV antigens and/or antibodies against FCoV are present in sample, a purple test line would appear in the result window.

The highly selective recombinant FCoV antigens and antibodies to FCoV are used as a capture and detector in the assay. These antigens and antibodies are capable of detecting FCoV antibodies and antigens in feline samples with high accuracy.

- Materials provided

| Odczynnik                            | 5 tests/kit |
|--------------------------------------|-------------|
| Test VetExpert Rapid FCoV Ag/Ab      | 5           |
| Assay diluents tube                  | 5           |
| Disposable dropper                   | 5           |
| Disposable swab                      | 5           |
| FCoV Ab Assay diluents bottle (2 ml) | 1           |
| Disposable capillary tube            | 5           |
| Anticoagulant tube                   | 5           |
| Instruction for use                  | 1           |

◆ A dark color score line on the capillary tube is the indicator line for 10 µl

● Materials required, but not provided

- 1) Timer

- Precautions

- 1) The kit is for feline use only. Do not apply to other animals.
- 2) The test device is sensitive to humidity as well as heat. Perform the test immediately after removing the test device from the foil pouch.
- 3) Do not reuse test components.
- 4) Apply the assay diluents and sample vertically.
- 5) Do not touch the membrane in the result window of test device.
- 6) Do not use the test kit beyond the stated expiration date marked on the label.
- 7) Do not use the test kit if the pouch is damaged or the seal is broken.
- 8) Do not mix components from different lot numbers because the components in this kit have been quality control tested as standard batch unit.
- 9) All sample should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- 10) Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.

- Storage and Stability

- 1) Store the test kit at 2–30°C. DO NOT FREEZE.
- 2) Do not store the test kit in direct sunlight.
- 3) The test kit is stable within the expiration date that marked on the package label.

- Collection and Preparation of Sample

#### [FCoV Ag]


- 1) Feline feces swab should be used for this test.

2) The samples should be tested immediately after collection.

3) If samples are not tested immediately, they should be stored at 2–8°C for 24 hours.

For longer storage, freeze at -20°C or below. Frozen samples should be brought to room temperature (15–30°C) prior to use.

4) The amount of fecal swab may affect the results. It is required to follow the swab amount of feces as shown in the picture on the left side.

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