VetExpert Rapid Test FIV Ab/FeLV Ag

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One-step Feline Immunodeficiency virus Antibody and Feline Leukemia virus Antigen Test

VetExpert Rapid FIV Ab/FeLV Ag Test

• Principles

The VetExpert Rapid FIV Ab/FeLV Ag Test is a chromatographic immunoassay for the qualitative detection of Feline Leukemia virus antigen and Feline Immunodeficiency virus antibody in feline serum, plasma or whole blood. The VetExpert Rapid Feline Leukemia virus Ag/Feline Immunodeficiency virus Ab Test Kit has the letters "T" and "C" as the Test line and Control line on the surface of the device. Both the test line and control line in the result window are not visible before applying any samples. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working. A purple test line will be visible in the result window if there is enough Feline Leukemia virus antigen and/or Feline Immunodeficiency virus antibody in the specimen.

The specially selected Feline Immunodeficiency virus antigen and Feline Leukemia virus antibody are used in the test both capture and detector materials. These enable the VetExpert Rapid FIV Ab/FeLV Ag Test to identify Feline Leukemia virus antigen and Feline Immunodeficiency virus antibody in feline serum, plasma or whole blood with a high degree of accuracy.

• Materials provided (2, 5, 10 test/kit)

1) 2, 5, 10 VetExpert Rapid FIV Ab/FeLV Ag Test

2) 1 or 2 assay diluents (3 ml)

3) 2, 5, 10 disposable capillary tubes

4) 2, 5, 10 anticoagulant bottles

5) Introduction for use

◆ A dark color score line on the capillary tube is the indicator line for 10 µl

• Precautions

1) For veterinary diagnostic use only.

2) For best results, strict adherence to the instructions is required.

3) All specimens should be handled as being potentially infectious.

4) Do not open or remove the test kits from their individually sealed pouches until immediately after their use.

5) Do not use the test kit if the pouch is damaged or the seal is broken.

6) Do not reuse test kits.

7) All reagents must be at room temperature before running the assay.

8) Do not use reagents beyond the stated expiration date marked on the label.

9) The components in this kit have been quality control tested as a standard batch unit. Do not mix components from different lot numbers.

• Storage and Stability

The kit can be stored at room temperature (2-30°C) or refrigerated. The test kit is stable through the expiration date marked on the package label. DO NOT FREEZE. Do not store the test kit in direct sunlight.

• Specimen Collection and Preparation

1) The test should be performed using serum, plasma, or whole blood.

2) [Whole blood] Collect on anticoagulated blood sample in EDTA, heparin or citrate using standard clinical laboratory procedures. Anticoagulated whole blood samples should be tested within 24 hours of drawing. If delays are expected, samples should be stored either on ice or refrigerated (2-7°C), but should not be frozen. If anticoagulated whole blood samples cannot be tested within this period of time, separate plasma by centrifugation and store as described in the next section.

• In-vitro-Diagnostikum

1) Test the specimen using serum, plasma, or whole blood.

2) [Whole blood] Collect on anticoagulated blood sample in EDTA, heparin or citrate using standard clinical laboratory procedures. Anticoagulated whole blood samples should be tested within 24 hours of drawing. If delays are expected, samples should be stored either on ice or refrigerated (2-7°C), but should not be frozen. If anticoagulated whole blood samples cannot be tested within this period of time, separate plasma by centrifugation and store as described in the next section.

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