
# VETEXPERT BG Vet

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# VETEXPERT BG Vet

## BLOOD GLUCOSE MONITORING SYSTEM FOR DOGS AND CATS.

For veterinary use only.

### IMPORTANT SAFETY PRECAUTIONS

#### Read Before Use

VETEXPERT BG Vet Blood Glucose Monitoring System consists of three main products: the meter, test strips and control solutions. These products have been designed, tested, and proven to work together as a system to provide accurate blood glucose test results.

**All products are sold separately.** Only use VETEXPERT BG Vet Blood Glucose Test Strips and control solutions with VETEXPERT BG Vet Blood Glucose Monitoring System.

#### Intended Use

This system is intended to be used for the quantitative measurement of glucose (sugar) in fresh whole blood samples. It should not be used as a single test for diagnosis or screening of diabetes mellitus. If anticoagulant is needed, please use heparin. We do NOT recommend using EDTA or NaF. System allows animal's blood glucose levels to be measured by the pet's owners at home and by veterinarian professionals in clinical settings as an aid to monitoring the effectiveness of diabetes control in cats and dogs. This system provides you with plasma equivalent results and are displayed in milligrams of glucose per deciliter of blood (mg/dL) or in millimoles of glucose per liter of blood (mmol/L).

#### APPEARANCE AND KEY FUNCTIONS OF THE METER

1. Test Strip Slot
2. Test Strip Ejector
3. Display Screen
4. M Button (M)
5. Data Port
6. SET Button (S)
7. Battery Compartment

#### NOTE

The meter will turn off automatically after 180 seconds without any action or you can press and hold the M for 3 seconds to turn off the meter

#### METER DISPLAY

1. Low Battery Symbol
2. Test Strip Symbol
3. Blood Drop Symbol
4. Code
5. Buzzer
6. Test Result
7. Day Average
8. Alarm Symbol
9. Date
10. Time
11. Measuring Mode
12. Control Solution Mode
13. Memory Mode
14. Error Message
15. Measurement Unit


#### VETEXPERT BG Vet STRIPS

Your system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole on the test strip and is automatically drawn into the reaction cell where the reaction takes place.

#### SETTING THE METER AND DELETING THE MEMORY


Start with the meter off (no test strip inserted). Press S.

##### STEP 1. Set the date


##### STEP 2. Setting the time format.

Press M to select the desired time format – 12h or 24h. Press S.

##### STEP 3. Setting the time

With the hour flashing, press M until the correct hour appears. With the minute flashing, press M until the correct minute appears. Press S.

##### STEP 4. Setting the unit of measurement

Press M to switch between mg/dL and mmol/L. Press S.


##### STEP 5. Setting the buzzer

With the buzzer displays, press M to switch between "On" and "OFF". Press S.

##### STEP 6. Deleting the memory

With "dEL" and "dE" on the display, press M and select "no" to keep the results in memory then press S to skip. To delete all the results, press M and select "yes" to delete all the memory records.

##### STEP 7. Setting the reminder alarm

Your meter has four reminder alarms. The meter will display "OFF" and "On".

If you don't want to set an alarm, press S to skip this step. Or press M to select "On", then press S.

With the hour/minute flashing in sequence, press M to select the correct hour/minute. Press S and go to the next alarm setting.

#### THE TWO MEASURING MODES

The meter provides you with two modes for measuring, General and QC.

MODES	USE WHEN
Gen	any time of day without regard to time since last meal
QC	testing with the control solution

You can switch between each mode by:

1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display code no., Gen and a flashing "d".

2. Press M to switch to QC mode.

#### BEFORE TESTING

##### Calibration

You must calibrate the meter every time you begin to use a new vial of test strips by setting the meter with the correct code. Test results may be inaccurate if the code number displayed on the monitor does not match the number printed on the strip vial.

##### How to Code Your Meter

1. Insert the code strip when the monitor is off. Wait until the code number appears on the display.

**Note:** Make sure the code number on display, code strip, and test strip vial are the same.

2. Remove the code strip, the display will show "OK". This tells you that the meter has finished coding and is ready for blood glucose testing.

##### Checking the Code Number

You need to make sure that the code number displayed on the meter matches the number on the test strip vial before you proceed. If it matches, you can proceed with your test. If the codes do not match, please stop testing and repeat the calibration procedure. If the problem persists contact Customer Service for help.

**NOTICE:** The code number on the display is only for your reference, it may not be actual code for this meter.

##### WARNING

It is important to make sure that the LCD displayed code is the same as the code on the test strip vial before testing. Failure to do so will get inaccurate results. If the LCD displayed code is not the same as the code on the test strip vial, and the code number cannot be updated, please contact Customer Service for assistance.

#### CHECKING THE SYSTEM WITH CONTROL SOLUTIONS

Our Control Solution contains a known amount of glucose that reacts with test strips. By comparing the result of your control solution test with the expected range printed on the test strip vial label, you can check that the meter and the test strips are working together as a system and that you are performing the test correctly. It is very important that you perform this simple check routinely to make sure you get accurate results.

##### How to perform a control solution test

Insert a test strip into the test slot with the contact bars end first and facing up. (Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.) The meter turns on automatically and displays code no., Gen and a flashing "d".

When the "QC" symbol appears on the display, press M and "QC" will appear on the display. When the "QC" sign is displayed, the meter will store your test result in the memory under QC. If you decide not to perform a control solution test, press M again and the "QC" sign will disappear.

##### CAUTION


Every time you perform a control solution test you must enter "QC" test mode so that the test result will NOT mix with the blood glucose TEST RESULTS stored in the meter memory. Failure to do so will mix up the blood glucose test results with the control solution test results in memory.

##### STEP 2. Apply control solution

Shake the control solution vial thoroughly before use. Squeeze out the first drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto a strip.

##### STEP 3. Read and compare the results

When the "QC" symbol appears on the display, press M and "QC" will appear on the display. When the "QC" sign is displayed, the meter will store your test result in the memory under QC. If you decide not to perform a control solution test, press M again and the "QC" sign will disappear. After the meter counts to 0, the result of the control solution test will appear. Compare the result with range printed on the test strip vial. The result should fall within this range.

#### Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do NOT test your animal's blood glucose level. Contact the local customer service for help.

#### TESTING WITH BLOOD SAMPLE

##### Preparing the Lancing Device for Blood Testing (not included)

Please follow the instructions in the lancing device insert for collecting a blood sample.

##### STEP 1. Insert the test strip to turn on the meter


Wait for the meter to display code no., Gen and a flashing "d".

##### STEP 2. Select the appropriate measuring mode by pressing M

For selecting the measurement mode, please refer to the "TWO MEASURING MODES".

##### STEP 3. Apply blood sample

Use the pre-set lancing device or a sterile, single-use injection needle to prick the skin of dog's marginal ear or cat's paw pad. Obtain a drop of blood at least 1.0 µL of volume.

Gently apply the drop of blood to the absorbent hole of the test strip at a tilted angle. Confirmation window should be completely filled if enough blood sample has been applied. Do NOT remove the paw of dog or cat until you hear a beep sound.

If the confirmation window is not filled completely before the meter begins to count down, do not add more blood to the test strip. Discard the test strip and start again. If you have trouble filling the test strip, please contact your veterinarian or the local customer service for assistance.

#### NOTE


If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.

##### STEP 4. Get result

The result of your pet's blood glucose test will be stored in the memory automatically.

##### STEP 5. Display hematocrit

To display the hematocrit level, press M and the result appeared. The hematocrit level will be stored in the memory.

#### WARNING

1. Please do not change your treatment on pets based on the result without first consulting your veterinarian.
2. Turn the meter off by removing the test strip. Discard the used test strip and lancet carefully according to your local regulations.

#### METER MEMORY


Your meter stores the 400 most recent blood glucose test results along with respective HCT values, dates and times in its memory. You can review the test results with these easy steps.

##### Reviewing Test Results

##### STEP 1. Enter the memory mode


With the meter turned off, press M twice. The first test result will appear, indicating that you are in the memory mode.

When using the meter for the first time or if the results have been deleted, "..." will appear, indicating that there are no test results in the memory.


##### STEP 2. Recall test results

Press M again, the hematocrit level for the first reading will show on the display. If you continue to press M, you can then review the last 400 results in the memory.

##### STEP 3. Exit the memory mode


After displaying the last test result in memory, press the M again. The meter displays "End" and then turns off.

##### Reviewing Day Average Results


##### STEP 1. Enter day average mode

With the meter turned off, press the M once. When "..." appears on the display, press M for 2 seconds to view day average result.

##### STEP 2. Recall test results average

Your 7-day average result will appear first. Press M to review 14-, 21-, 28-, 60- and 90-day average results stored in memory.

##### STEP 3. Exit the memory mode

Press the M for three seconds to turn off the meter.

#### DOWNLOADING RESULTS ONTO A COMPUTER

##### Data Transmission via Cable (not included)

You can use the meter with an RS232 cable and the Health Care Software System to view test results on your personal computer. To learn more about the Health Care Software System or to obtain a RS232 cable separately, please contact the local customer services for assistance.

#### WARNING

While the meter is connecting to the PC, it will be unable to perform a blood glucose test.

#### BATTERY

Your meter comes with one 3V CR2032 lithium battery.

If the low battery symbol "b" appears on the screen indicating that the battery is low and it is time to change the battery.

The "b" appears with E-b. Error and low:

The power is not enough to do a test. Please change the battery immediately.

##### Replacing the Battery