VET EXPERT BG Vet

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VET EXPERT BG Vet

BLOOD GLUCOSE TEST STRIPS FOR DOGS AND CATS

Use only with VET EXPERT BG Vet / BG Vet Pro Blood Glucose Monitoring System to test for glucose in dogs and cats.
For veterinary use only.

WARNINGS

- For veterinary use only.
- For single use only.
- Please read this sheet and the VET EXPERT BG Vet / BG Vet Pro Operation Instruction before you use this test strip. Use only VET EXPERT BG Vet Test Strips with the VET EXPERT BG Vet / BG Vet Pro Blood Glucose Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Severe dehydration and excessive water loss may cause false low results. If severe dehydration is suspected, consult a veterinarian/veterinary surgeon immediately.
- Do not use the system to test an animal's blood until a control solution test result is obtained that is within the recommended range. The recommended control solution range is on the test strip vial.
- Keep test strips away from small children. If swallowed, consult a doctor immediately for advice.

INTENDED USE

VET EXPERT BG Vet test strips, when used together with the VET EXPERT BG Vet / BG Vet Pro blood glucose monitoring system, allow pet's blood glucose levels to be measured by the pet's owner at home and by veterinarian professionals in clinical settings as an aid to monitoring the effectiveness of diabetes control in cats and dogs. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh whole blood samples. It should not be used as a single test for diagnosis or screening of diabetes mellitus. Professionals may test with capillary and venous blood sample from paw pad and marginal ear.

LIMITATIONS

- Hematocrit: The hematocrit level is limited to between 0% and 70%.
- This test strip is used for testing fresh capillary and venous blood.
- Altitude Ects: Altitudes up to 10,742 feet (3,275 m) do not affect test results.
- Use only heparin for anticoagulation of whole blood.


STORAGE AND HANDLING

IMPORTANT

Do not use the test strips if they have expired.


- Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it.
- Store the test strips in a cool, dry place between 2°C and 30°C (35.6°F and 86°F) and below 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.
- Store the test strips in their original vial ONLY. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Use each test strip immediately after taking it out of the vial or individual foil packet. Close the vial immediately after taking out a strip.
- Keep the vial closed at all times.
- Do not bend, cut, or alter the test strip.

STRIP APPEARANCE

ATTENTION

Test results might be wrong if the contact bar is not fully inserted into the test slot. The front side of the test strip should face up when inserting the test strip.

TESTING ANIMAL'S BLOOD GLUCOSE

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTING.

STEP 1.

Insert the test strip fully into the slot of the meter until it will go no further. When the strip is fully inserted, the meter will do several self-checks.

Please make sure you have calibrated your meter with the code strip according to the instructions in the Operation Instruction before testing.

STEP 2.

Collect a blood sample for about 1.0 µL with the test strip. A sufficient quantity of blood is required for the test to provide accurate results. Touch the blood drop with the absorbent hole of the test strip, and wait until the confirmation window is fully covered. Do NOT apply a smeared blood sample. The meter will start counting down.

STEP 3.

After a few seconds, the meter will display your pet's blood glucose level. The last reading will be automatically saved in the meter. Turn it off by removing the test strip and throw away the used test strip.

Make sure you use the code strip and calibrate your meter before using for the first time or changing a new vial of test strip. Please refer to your Operation Instruction for more information.

The used test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

READING RESULTS

The blood glucose readings deliver plasma equivalent results and are displayed in milligrams of glucose per deciliter of blood (mg/dL) or in millimoles of glucose per liter of blood (mmol/L). The measurement range of this meter is 10-600 mg/dL (0.5-33.3 mmol/L). Please consult a veterinarian/veterinary surgeon for the target glucose values for each animal.

Questionable or inconsistent results

- If the test results are unusual or inconsistent:
 - Make sure the confirmation window of the test strip is completely filled with blood.
 - Check the expiration date of the test strips.
 - Check the performance of your meter and test strip with the control solutions.
 - Make sure your monitor has correct coding, and the code is the same as the code on the vial or individual foil packet you are using.

PLEASE NOTE

Low or high blood glucose readings can indicate a potentially serious medical condition. If the blood glucose reading is unusually low or high repeat the test with a new test strip. If the reading is not consistent with any visible symptoms, contact a veterinarian/veterinary surgeon and follow his or her treatment advice.

QUALITY CONTROL TESTING

Our control solutions contain a known amount of glucose that can react with test strips. You can check the performance of meter, test strip and your technique by comparing the control solution results with the range printed on the label of test strip vial. Checking regularly can ensure the test results are accurate. Please refer to the Operation Instruction for complete testing instructions.

IMPORTANT

The reference range of the control solutions may vary with each new vial or package of test strips. Make sure you check the range on the label of your current vial.

CHEMICAL COMPONENTS

| | |
|--|-----|
| Glucose dehydrogenase (<i>E. coli</i>) | 8% |
| Electron shuttle | 55% |
| Enzyme protector | 8% |
| Non-reactive ingredients | 29% |

ADDITIONAL INFORMATION FOR VETERINARY PROFESSIONALS

Veterinary professionals may use test strips to test capillary and venous whole blood. Use fresh whole blood samples only.

Sample Size: 1.0 L

Reaction Time: 5 seconds

System Measurement Range: 10 mg/dL to 600 mg/dL (0.5 mmol/L to 33.3 mmol/L)

Hematocrit Range: 0% to 70%

SYMBOL INFORMATION

| Symbol | Referent | Symbol | Referent |
|--------|---|--------|---------------------|
| | Consult instructions for use | | Manufacturer |
| | Caution, consult accompanying documents | | Do not reuse |
| | Use by/EXP | | Humidity limitation |
| | Temperature limitation | | RoHS Compliance |
| | Batch code | | CE mark |

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