# VET EXPERT BG Vet

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# VET EXPERT BG Vet

## BLOOD GLUCOSE CONTROL SOLUTION

Use only with Vet Expert Monitoring Systems.  
For veterinary use only.

### INTENDED USE

VET EXPERT Blood Glucose Control Solution is intended for use with VET EXPERT Blood Glucose Strips. It is used to check that the meter and test strips are working together properly. It is important to do control solution tests routinely to ensure getting accurate results.

### SUMMARY AND EXPLANATION

VET EXPERT Blood Glucose Control Solution contains a known amount of glucose (NORMAL and HIGH) that reacts with test strips.

Operation of the control solution test is similar to blood glucose test except that you use blood glucose control solution instead of a drop of blood. The control solutions should provide results within an expected range printed on the test strip vial label.

#### IMPORTANT

The control range may change with each new vial of test strips. Always use the control range on the label of your current vial of test strips.

#### How often should the control solution test be performed?

- When you use this system to test the blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the expected range, you are ready to test the blood.
- For routinely check the meter and test strips, perform a single test for each level of control solution at least once a week.

#### When should the control solution test be performed?

- Whenever you suspect that the meter or test strips are not working properly.
- When blood glucose test results are inconsistent with clinical condition of the animal, or when results are inaccurate.
- When you want to practice the testing process.
- If you drop the meter or think you may have damaged the meter.

#### CAUTION

The control solution range printed on the test strip vial is for control solution only. It is used to test meter and test strip performance. It is not a recommended range for blood glucose level.

### CHEMISTRY

Blood Glucose Control Solution is a red solution that contains less than 0.30% D-glucose. The solution does not contain any human-derived or biological materials.

Ingredient	Percent
Active ingredient: D-glucose	0.03-0.30%
Inert ingredients	>99.70%

### STORAGE AND HANDLING

- For veterinary use only.
- For in vitro diagnostic use only.
- Use only VET EXPERT Blood Glucose Control Solution with VET EXPERT Glucose Test Strip.
- Check the expiration date on the vial. Do not use if expired.
- Tightly close the control solution and store at room temperature between 2°C and 30°C. Do not freeze.
- Use only for 90 days after first opening. Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.

### TEST PROCEDURE

To do a test with VET EXPERT Blood Glucose Control Solution, follow the Operation Instruction that came with the system.

Accurate results depend upon following instructions carefully. Follow below technical pointers whenever you test with VET EXPERT Blood Glucose Control Solution:

1. It is recommended that the meter, control solution, and test strips come to room temperature (20°C-25°C) before testing.
2. Shake the control solution well before using.
3. Discard the first drop of control solution you squeeze out.
4. Wipe off the dispenser tip to keep it clean.

### EXPECTED RESULTS

At room temperature, VET EXPERT Blood Glucose Control Solution test results should fall within the expected range printed on the test strip vial label.

If test results fall outside the expected range printed on the test strip vial, please repeat the test using a new strip. Results that fall outside the expected range may indicate:

- An error in the way you are doing the test.
- Expired test strips and control solutions.
- A problem with the meter.
- Control solution that is too warm or too cool.
- The meter code does not match the code printed on the test strip vial.
- Failure to shake the control solution vial vigorously.

If the repeat result is still outside the range, the system may not be working properly. Do not use the system to test your blood. Please contact the dealer for service.

### SYMBOL INFORMATION

Symbol	Referent	Symbol	Referent
	In vitro diagnostic medical device		Batch code
	Consult instructions for use		Manufacturer
	Temperature limitation		CE mark
	Use by /EXP		Caution, consult accompanying documents

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